# Sippi® Urinometro Digitale

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### A new automatic urinometer is more exact and higher rated than a manual urinometer

Anton Eklund, Martin Slettengren, Jan van der Linden  
Department of Cardiothoracic Surgery and Anesthesiology, Karolinska University Hospital, Stockholm, Sweden

#### Conclusions

The automatic urinometer has a lower mean bias compared with the manual urinometer.

Possible imprecision due to temporal deviations using the manual urinometer is avoided by the continuous measurement of the automatic urinometer.

The automatic urinometer was rated higher in staff satisfaction, reflecting perception of higher reliability, ease of use and less time usage.

The features of the automatic urinometer may also indicate a favorable clinical impact in the normal ward, where hourly measurements with a manual urinometer usually are too time-consuming.

Karolinska Institutet.

Publicato su Critical Care, 2015.

Karolinska Institute  
Dept. Cardiothoracic Surgery  
and Anesthesiology  
Karolinska University Hospital  
Stockholm, Sweden

E-mail: antoneklund@hotmail.com  
martin.slettengren@karolinska.se  
jan.vanderlinden@ki.se  
Phone: +46 8 51770121

### Liquid-Infused Silicone As a Biofouling-Free Medical Material

Noah MacCallum,<sup>†</sup> Caitlin Howell,<sup>†,‡</sup> Philseok Kim,<sup>†</sup> Derek Sun,<sup>†</sup> Ronn Friedlander,<sup>‡</sup> Jonathan Ranisau,<sup>†</sup> Onye Ahanotu,<sup>†</sup> Jennifer J. Lin,<sup>†</sup> Alex Vena,<sup>†</sup> Benjamin Hatton,<sup>†</sup> Tak-Sing Wong,<sup>†</sup> and Joanna Aizenberg<sup>\*,†,‡,§</sup>

<sup>†</sup>Wyss Institute for Biologically Inspired Engineering, <sup>‡</sup>School of Engineering and Applied Sciences, and <sup>§</sup>Department of Chemistry and Chemical Biology, Harvard University, Cambridge, Massachusetts, United States

#### Supporting Information

**ABSTRACT:** There is a dire need for infection prevention strategies that do not require the use of antibiotics, which exacerbate the rise of multi- and pan-drug resistant infectious organisms. An important target in this area is the bacterial attachment and subsequent biofilm formation on medical devices (e.g., catheters). Here we describe nonfouling, lubricant-infused slippery polymers as proof-of-concept medical materials that are based on oil-infused polydimethylsiloxane (iPDMS). Planar and tubular geometry silicone substrates can be infused with nontoxic silicone oil to create a stable, extremely slippery interface that exhibits exceptionally low bacterial adhesion and prevents biofilm formation. Analysis of a flow culture of *Pseudomonas aeruginosa* through untreated PDMS and iPDMS tubing shows at least an order of magnitude reduction of biofilm formation on iPDMS, and almost complete absence of biofilm on iPDMS after a gentle water rinse. The iPDMS materials can be applied as a coating on other polymers or prepared by simply immersing silicone tubing in silicone oil, and are compatible with traditional sterilization methods. As a demonstration, we show the preparation of silicone-coated polyurethane catheters and significant reduction of *Escherichia coli* and *Staphylococcus epidermidis* biofilm formation on the catheter surface. This work represents an important first step toward a simple and effective means of preventing bacterial adhesion on a wide range of materials used for medical devices.

**KEYWORDS:** nonfouling material, biofilm prevention, slippery surfaces, nosocomial infection, medical materials

Bacteria-Free Treated Side | Contaminated Untreated Side

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SIM Italia Srl